



# OIT CONNECT

A GFT COMPANY



Designed by OIT-Allergists, for OIT-Practices & OIT-Patients



# Oral Immunotherapy Software by the Leaders in OIT

Helping patients and physicians  
connect to enhance the OIT  
experience.



[www.oitconnect.com](http://www.oitconnect.com)



# Designed by OIT-Allergists for OIT-Practices & Patients



## Increase Practice Efficiency

Built-in protocols and workflows designed specifically for OIT by experienced leaders in the OIT field ensures you have everything you need for the most efficient OIT documentation.



## Enhance Communication

Keep patients engaged. With **two-way secure messaging** between providers, practice, and the connected patient application, patients can better manage their OIT treatment plan, document maintenance dosing, and communicate with the practice.



## Improve Quality of Life

The goal of OIT is to improve the quality of life of patients.

With the OITConnect™ **cloud-based platform**, our clients can access the system from anywhere to better manage OIT treatment, see results, and improve their own quality of life with a more efficient workflow.




- **Our mission is to help you help your patients**
- **Having OIT food products or protocols is only the beginning**
- **Get the right technology in OITConnect™ to help drive results and enhance safety and efficiency**



# What is OITConnect™

- **Easy, cloud-based, HIPAA-compliant software**
  - **Stand alone cloud-based**
  - **Integrate with your existing EHR**
  - **ModuleMD – One Click Access**
- **Enhance safety and efficiency of an OIT clinical practice**
- **Simple way to document everything with Food OIT**
- **Practice Portal**
- **Integrated Patient Portal**

# LOGIN – Easy, Cloud-based, HIPAA Compliant



The screenshot displays the OIT Connect website. The header includes the OIT Connect logo (A GFT COMPANY) and navigation links: [ Home ], About, Patients, Providers, and Contact. A 'LOGIN' button is highlighted in the top right corner with a yellow box and an arrow pointing to a larger, detailed view of the login form on the right. The main content area features a large image of a family (a man, a woman, and two children) sitting at a dining table with food. Overlaid on this image is the text: 'Better Quality of Life with Oral Immunotherapy (OIT)' and 'Food should be a source of nourishment and connection, not stress. OIT makes that possible.' Below this text is a green 'Learn More' button. The detailed login form on the right includes the OIT Connect logo, the text 'Sign in with your email address', an 'Email Address' input field, a 'Password' input field, a 'Forgot your password?' link, a blue 'Sign in' button, and a 'Don't have an account? Sign up now' link.

Registered users without an EMR integration can login through the OITConnect website at [www.oitconnect.com](http://www.oitconnect.com)

# PRACTICE SETUP



# FOODS – Preloaded Foods and Protocols – Single or Multi-Food OIT

We have set up the primary OIT Food categories, so you don't have to.

Peanut  
Almond  
Eggs  
Cashew  
Cow's Milk  
Wheat  
Walnut  
Soy  
Sesame Seeds  
Sesame Tahini

We can add  
as many foods and  
protocols  
as you need!

**OIT CONNECT** 15:37:04  
Medical Clinic Inc.  
DANDAMUDI SESHAGIRI, Detroit (EST)

**My Activities**  
Patient Chart Messages Reception Schedule Schedule Reports Finalize Notes Practice Setup

Select Patient Search Options

**Dashboard Practice Setup**

General Billing Schedule EMR Privileges Register IP Address OIT

**Setup OIT**

**Foods Food Types Solutions Protocols Reactions Medications Units Concentrations Manufacturers**


Search: Status: Active New

Name	Active
Peanut	Yes
Almond	Yes
Eggs	Yes
Cashew	Yes
Cow's Milk	Yes
Wheat	Yes
Walnut	Yes
Hazelnut	Yes
Soy	Yes
Sesame Seeds	Yes
Sesame Tahini	Yes

**Update Food**

Name: Peanut

Active: ☒







# PRELOADED OIT FOOD PRODUCTS

<b>OITCONNECT</b> 15:41:43 Medical Clinic Inc. DANDAMUDI SESHAGIRI, Detroit (EST)			
<b>My Activities</b> Patient Chart Messages Reception Schedule Schedule Reports Finalize Notes Practice Setup			
Select Patient Search Options			
<b>Dashboard Practice Setup</b>			
General Billing Schedule EMR Privileges Register IP Address OIT			
<b>Setup OIT</b>			
Foods Food Types Solutions Protocols Reactions Medications Units Concentrations Manufacturers			
Search: Status: Active New Food Type Inventory			
Name	Food	Type	Active
Distilled Water		Diluent	Yes
Kool Aid		Additive	Yes
Palforzia Peanut 240 mg capsule		Solute	Yes
Apple Juice		Diluent	Yes
Almonds	Almond	(1.00 mg) Protein	Yes
Almond Milk	Almond	Solute	Yes
Cashew Fragments	Cashew	Solute	Yes
Cashew	Cashew	Solute	Yes
Cashew Flour	Cashew	Solute	Yes
Chickpea flour	Chickpea	Solute	Yes
Chickpea whole	Chickpea	(534.00 mg) Protein	Yes
Coconut Flour	Coconut	(2.00 mg) Protein	Yes
Coconut Shredded.	Coconut	( mg) Protein	No
Coconut Shredded	Coconut	( mg) Protein	No
CoconutShredded	Coconut	( mg) Protein	No
Coconut Shredded	Coconut	( mg) Protein	Yes
Baked Milk Muffins	Cow's Milk	Solute	Yes
Horizon Organic Vit D Whole Milk	Cow's Milk	Solute	Yes
All Whites 100% Liquid Egg Whites	Eggs	Solute	Yes
Egglands' Best Liquid Egg Whites	Eggs	Solute	Yes
Baked Egg	Eggs	Solute	Yes
Egg whole	Eggs	(4995.00 mg) Protein	Yes
Hazelnut Fragments	Hazelnut	Solute	Yes
Hazelnut	Hazelnut	Solute	Yes
Hazelnut Flour	Hazelnut	Solute	Yes
Macadamia Fragments	Macadamia	Diluent	Yes
Horizon Organic Milk	Milk	Solute	Yes
0.5 mg Capsule	Palforzia Peanut	(0.50 mg) Protein	No
1 mg Capsule	Palforzia Peanut	(1.00 mg) Protein	No
Palforzia Peanut 1 mg Capsule	Palforzia Peanut	(1.00 mg) Protein	Yes
Palforzia Peanut 3 mg Capsule	Palforzia Peanut	(3.00 mg) Protein	Yes
Palforzia Peanut 6 mg Capsule	Palforzia Peanut	(6.00 mg) Protein	No
12 mg capsule	Palforzia Peanut	(12.00 mg) Protein	No

Distilled Water  
Kool Aid  
Almond Milk  
Almonds  
Cashew  
Cashew Flour  
Cashew Fragments  
Chickpea flour  
Chickpea whole  
Coconut Flour  
Coconut Shredded  
Coconut Shredded  
Coconut Shredded.  
Coconut Shredded  
Baked Milk Muffins  
Horizon Organic Vit D Whole Milk  
All Whites 100% Liquid Egg Whites  
Baked Egg  
Egg whole  
Egglands' Best Liquid Egg Whites  
Hazelnut  
Hazelnut Flour  
Hazelnut Fragments  
Macadamia Fragments  
Horizon Organic Milk  
Palforzia Peanut 240 mg capsule  
Palforzia Peanut 0.5 mg Capsule  
Palforzia Peanut 1 mg Capsule  
Palforzia Peanut 12 mg capsule  
Palforzia Peanut 120 mg  
Palforzia Peanut 160 mg  
Palforzia Peanut 20 mg capsule  
Palforzia Peanut 200 mg  
Palforzia Peanut 240 mg  
Palforzia Peanut 300 mg  
Palforzia Peanut 40 mg capsule  
Palforzia Peanut 80 mg

Palforzia Peanut 0.5 mg Capsule  
Palforzia Peanut 1 mg Capsule  
Palforzia Peanut 12 mg capsule  
Palforzia Peanut 120 mg capsule  
Palforzia Peanut 160 mg capsule  
Palforzia Peanut 20 mg capsule  
Palforzia Peanut 200 mg capsule  
Palforzia Peanut 240 mg capsule  
Palforzia Peanut 3 mg Capsule  
Palforzia Peanut 300 mg capsule  
Palforzia Peanut 40 mg capsule  
Palforzia Peanut 6 mg Capsule  
Palforzia Peanut 6 mg capsule  
Palforzia Peanut 80 mg capsule  
Palforzia Peanut 9 mg capsule  
1000 mg Peanut  
Palforzia Peanut 1.5 mg  
Peanut Flour  
Peanut Fragments  
Peanut Solution  
Whole Peanut  
Rye Bread  
RYE Flour  
Sesame Flour  
Sesame Seeds (Mortar & Pestle Ground)  
Sesame Seeds (Unhulled)  
Sesame (Unhulled)  
Sesame Tahini  
Soy Milk (Silk Organic)  
Test Food Type  
Walnut Fragments  
Walnut Meal (Flour)  
Walnuts  
Wheat Bread (4g/slice)  
Wheat Flour (Whole Wheat)



# INVENTORY – MANAGE EXPIRATION DATES AND STOCK

Manage your in-office Food Inventory for an easier mixing process



## Inventory Management

Manage OIT inventory for your practice with the OITConnect™ inventory management feature. Our unique OIT specific inventory solution Works with **FDA approved OIT products** and allows you to track manufacturers, formula, expiration dates, additives, and more.

The screenshot displays the OITConnect Medical Clinic Inc., Detroit (EST) interface. The 'Update Food Type Inventory' dialog box is open, showing the following fields:

- Lot#: 12345
- Manufacturer: America's Peanut Farm
- Quantity: 455.00 ml
- Effective Date: 10/25/2020
- Expiry Date: 10/31/2020

The background shows the main application window with a table of food types. The table has columns for Lot#, Manufacturer, Quantity, Effective Date, and Expiry Date. The first row shows Lot# 12345, Manufacturer America's Peanut Farmers, Quantity 455.00 ml, Effective Date 10/25/2020, and Expiry Date 10/31/2020.

# PRELOADED AND CUSTOMIZABLE SOLUTIONS

A full list of OIT solution set have already been added for easy OIT Ordering.

OITCONNECT Medical Clinic Inc., Detroit (EST)

Select Patient  Search Options

Dashboard Practice Setup X

General Billing Schedule EMR Privileges Register IP Address OIT


Setup OIT X

Foods Food Types Solutions Protocols Reactions Medications Units Concentrations Manufacturers

Search:  Status: Active  New Solution Inventory

Name	Concentration	Volume (ml)
Almond (4)		
Almond Solution A	1 mcg/ml	100.00
Almond Solution B	10 mcg/ml	100.00
Almond Solution C	100 mcg/ml	100.00
Almond Solution D (Stock)	1.04 mg/ml	100.00
Cashew (4)		
Cashew Solution A	5 mcg/ml	100.00
Cashew Solution B	50 mcg/ml	100.00
Cashew Solution C	500 mcg/ml	100.00
Cashew Solution D (Stock)	5 mg/ml	250.00
Cow's Milk (3)		
Cow's Milk Soln A	100 mcg/ml	100.00
Cow's Milk Soln B	1 mg/ml	200.00
Cow's Whole Milk Horizon Organic	Whole Milk Stock	100.00
Eggs (3)		
Liq Egg White (Stock)	111 mg/ml	10.00
Liq Egg White Solution A	100 mcg/ml	100.00
Liq Egg White Solution B	1 mg/ml	111.00
Hazelnut (4)		
Hazelnut Solution A	15 mcg/ml	40.00
Hazelnut Solution B	150 mcg/ml	40.00
Hazelnut Solution C	1.5 mg/ml	40.00
Hazelnut Solution D (Stock)	15 mg/ml	100.00
Peanut (7)		
PN Solution AA	0.25 mcg/ml	40.00
PN Solution A	2.5 mcg/ml	40.00
PN Solution B	25 mcg/ml	40.00
PN Solution C	250 mcg/ml	40.00
PN Solution D	2.5 mg/ml	40.00
PN Solution E	25 mg/ml	40.00
PN Solution F (Stock)	250 mg/ml	50.00

Name:

Food:  

Concentration:

Volume:  ml


Protein:  mg

Start:  Active:

Solution Details - Peanut

Solutions / Food Types	Quantity	Unit	Protein (mg)	Add
Peanut Flour	12.250	grams		
Distilled Water	50.000	ml		<input type="button" value="+"/>

Name:

Food:  

Concentration:

Volume:  ml

Protein:  mg

Start:  Active:

Solution Details - Peanut

Solutions / Food Types	Quantity	Unit	Protein (mg)	Add
PN Solution B 25 mcg/ml	4.000	ml		
Kool Aid	5.000	grams		
Distilled Water	36.000	ml		<input type="button" value="+"/>

Manage all your OIT orders. You can **Order, Mix, Confirm**, document **Medications**, and **Administer** OIT doses from this screen.

**OITCONNECT** 19:48:34 **My Activities**  
 Medical Clinic Inc.  
 DANDAMUDI SESHAGIRI, Detroit (EST)





Patient Chart Messages Reception Schedule Schedule Reports Finalize Notes Practice Setup

abraham, sarah Search Options Patient: abraham, sarah Account #: 97800 DOB: 11/10/2017 Age: 3 yrs, 7 mo Sex: F Logout

Dashboard

Daily Patient Triage Office Communication Messages Patient Portal **OIT Orders**

New Edit Cancel Order Confirm Order Dose Administer Mix Medication Finalize View

Order Id	Order Date	Appointment Date	Order Status	Account #	Patient Name	Date Of Birth	Food	Image	Provider	Location	Confirmed By	Confirmed Date	Finalized By	Finalized Date
145	06/03/2021	06/07/2021	Dose Administered	97757	OITConnect, Patient	07/13/1972	Peanut		Giri Dandamudi, M.D	Detroit	Dandamudi, Seshagiri	06/03/2021 11:40 AM		
144	06/03/2021	06/03/2021	Dose Administered	97739	Demo, Mark	01/01/1998	Peanut		Giri Dandamudi, M.D	Detroit	Dandamudi, Seshagiri	06/03/2021 11:23 AM		
142	05/05/2021	05/05/2021	Ordered	97800	abraham, sarah	11/10/2017	Cashew		Douglas Jones	Detroit				
139	05/05/2021	05/05/2021	Ordered	97800	abraham, sarah	11/10/2017	Cashew		Douglas Jones	Detroit				


To order OIT for a patient:

**Patient Chart > OIT > Orders**

Select the **New**.

Here you will need to select the Provider, Food, and Protocol from the dropdown menu.

Review the auto populated patient treatment plan. When complete, select **Save**.



Medical Clinic Inc.  
DANDAMUDI SESHAGIRI, Detroit (EST)

20:01:11

My Activities

Patient Chart Messages Reception Schedule Schedule Reports Finalize Notes Practice Setup

OITConnect, Patient Test Search Options Patient: OITConnect, Patient Test Account #: 97757 DOB: 07/13/1972 Age: 48 yrs, 11 mo Sex: M More

Park Tabs ? 1 Logout


Dashboard Patient Chart X

Demographics View Notes Visit Notes Messages OIT

Patient Dashboard Orders X

New Edit Cancel Order Confirm Order Dose Administer Mix Medication Finalize

Order Id	Order Date	Appointment Date	Order Status	Account #	Patient Name	Date Of Birth	Food	Image	Provider	Location	Confirmed By	Confirmed Date	Finalized By	Finalized Date
145	06/03/2021	06/07/2021	Dose Administered	97757	OITConnect, Patient	07/13/1972	Peanut		Giri Dandamudi, M.D	Detroit	Dandamudi, Seshagiri	06/03/2021 11:40 AM		



Medical Clinic Inc.  
DANDAMUDI SESHAGIRI, Detroit (EST)

20:02:15

My Activities

Patient Chart Messages Reception Schedule Schedule Reports Finalize Notes Practice Setup

OITConnect, Patient Test Search Options Patient: OITConnect, Patient Test Account #: 97757 DOB: 07/13/1972 Age: 48 yrs, 11 mo Sex: M More

Park Tabs ? 1 Logout

Dashboard Patient Chart X

Demographics View Notes Visit Notes Messages OIT

Patient Dashboard Orders X

Add Order

Order Date: 06/15/2021 Provider: ATUL SHAH Location: Detroit

Food: Peanut Protocol: Peanut

Seq#	Dosing	Location	Solutions / Food Types	Dose	Unit	Protein (mg)	Frequency	Duration(In Days)	Add	Delete
1	DAY ONE	In Office	PN Solution A 2.5 mcg/ml	2.00	ml		every 20 minutes		+	-
2	DAY ONE	In Office	PN Solution A 2.5 mcg/ml	4.00	ml		every 20 minutes		+	-
3	DAY ONE	In Office	PN Solution B 25 mcg/ml	1.00	ml		every 20 minutes		+	-
4	DAY ONE	In Office	PN Solution B 25 mcg/ml	2.00	ml		every 20 minutes		+	-
5	DAY ONE	In Office	PN Solution B 25 mcg/ml	4.00	ml		every 20 minutes		+	-
6	DAY ONE	In Office	PN Solution C 250 mcg/ml	1.00	ml		every 20 minutes		+	-
7	DAY ONE	In Office	PN Solution C 250 mcg/ml	2.00	ml		every 20 minutes		+	-
8	DAY ONE	In Office	PN Solution C 250 mcg/ml	4.00	ml		every 20 minutes		+	-
9	HOME DOSE	At Home	PN Solution C 250 mcg/ml	4.00	ml		once a day	6	+	-
10	UPDOSE	In Office	PN Solution D 2.5 mg/ml	1.00	ml		once a day	1	+	-
11	HOME DOSE	At Home	PN Solution D 2.5 mg/ml	1.00	ml		once a day	6	+	-
12	UPDOSE	In Office	PN Solution D 2.5 mg/ml	2.00	ml		once a day	1	+	-
13	HOME DOSE	At Home	PN Solution D 2.5 mg/ml	2.00	ml		once a day	6	+	-
14	UPDOSE	In Office	PN Solution D 2.5 mg/ml	4.00	ml		once a day	1	+	-



Customized Treatment Plans Through the OITConnect™ food immunotherapy ordering system, practices can customize a treatment plan schedules that are appropriate for each patient's OIT journey. Review printable PDF reports of clinical notes and completed OIT protocols for trending throughout each patient's treatment.



# MIXING: OIT ORDERS

## Dashboard > OIT Orders

Select the **checkbox** next to the orders you would like to mix. All the orders selected should be for the same **Food**.

Select **Mix**.

**OITCONNECT** 20:10:53

Medical Clinic Inc.  
DANDAMUDI SESHAGIRI, Detroit (EST)

**My Activities**

Patient Chart Messages Reception Schedule Schedule Reports Finalize Notes Practice Setup

Demo, Mark Search Options Patient: Demo, Mark Account #: 97739 DOB: 01/01/1998 Age: 23 yrs, 5 mo Sex: M More.

Dashboard

Daily Patient Triage Office Communication Messages Patient Portal **OIT Orders**

New Edit Cancel Order Confirm Order Dose Administer Mix Medication Finalize

	Order Id	Order Date	Appointment Date	Order Status	Account #	Patient Name	Date Of Birth	Food	Image	Provider	Location	Confirmed By	Confirmed Date
<input checked="" type="checkbox"/>	144	06/03/2021	06/03/2021	Dose Administered	97739	Demo, Mark	01/01/1998	Peanut		Giri Dandamudi, M.D	Detroit	Dandamudi, Seshagiri	06/03/2021 11:23 AM
<input checked="" type="checkbox"/>	146	06/15/2021	06/07/2021	Ordered	97757	OITConnect, Patient	07/13/1972	Peanut		Giri Dandamudi, M.D	Detroit		
<input checked="" type="checkbox"/>	145	06/03/2021	06/07/2021	Dose Administered	97757	OITConnect, Patient	07/13/1972	Peanut		Giri Dandamudi, M.D	Detroit	Dandamudi, Seshagiri	06/03/2021 11:40 AM
<input type="checkbox"/>	142	05/05/2021	05/05/2021	Ordered	97800	abraham, sarah	11/10/2017	Cashew		Douglas Jones	Detroit		
<input type="checkbox"/>	139	05/05/2021	05/05/2021	Ordered	97800	abraham, sarah	11/10/2017	Cashew		Douglas Jones	Detroit		



When completed, select **Submit**.

[www.oitconnect.com](http://www.oitconnect.com)



# SOLUTION MIXING – INVENTORY – LABELS – PATIENT LABELS

Manage your in-office Solution inventory including Lot #, Volume, Expiry Date, and print bottle labels.

The screenshot displays the OITConnect software interface. The top navigation bar includes 'My Activities' with links for Patient Chart, Messages, Reception, Schedule, Schedule Reports, Finalize Notes, and Practice Setup. Below this is a 'Select Patient' dropdown and a 'Search Options' field. The main menu on the left includes 'Dashboard', 'Practice Setup', 'General', 'Billing', 'Schedule', 'EMR', 'Privileges', 'Register IP Address', and 'OIT'. The 'OIT' section is expanded, showing 'Setup' and 'OIT X' tabs. Under 'OIT X', there are sub-tabs for 'Foods', 'Food Types', 'Solutions', 'Protocols', 'Reactions', 'Medications', 'Units', 'Concentrations', and 'Manufacturers'. The 'Solutions' tab is active, showing a table for 'Solution Inventory - Almond Solution A'. The table has columns for Lot#, Volume, Effective Date, Expiry Date, and Active. A row is visible with Lot# 123, Volume 6.00 ml, Effective Date 03/24/2021, Expiry Date 04/24/2021, and Active Yes. To the right of the table is a 'Print Label' button. Below the table is a 'New' button. On the right side of the interface, there is a 'Update Solution Inventory' form. It includes fields for Lot# (123), Volume (6.00 ml), Effective Date (03/24/2021), Expiry Date (04/24/2021), and Active (checked). Below these fields is a 'Solution Details Almond' table with columns for Solutions / Food Types, Quantity, Mix Quantity, Lot#, and Expiry Date. The table has two rows: 'Almond Solution B 10 mc...' with Quantity 10.00 ml and Mix Quantity 0.60, and 'Distilled Water' with Quantity 90.00 ml and Mix Quantity 5.40. At the bottom right of the interface, there is a box containing the text: 'Almond Almond Solution A 1 mcg/ml Lot#123 Expireis:04/24/2021'.

Lot#	Volume	Effective Date	Expiry Date	Active
123	6.00 ml	03/24/2021	04/24/2021	Yes

Solutions / Food Types	Quantity	Mix Quantity	Lot#	Expiry Date
Almond Solution B 10 mc...	10.00 ml	0.60		
Distilled Water	90.00 ml	5.40		

Almond  
Almond Solution A 1 mcg/ml  
Lot#123 Expireis:04/24/2021



## OIT Mixing

With OITConnect™, practices can manage sequences and access information for mixing all in one place. Additionally, **content labels can be conveniently printed** to enable safer and more efficient mixing workflows.





# PRELOADED, YET CUSTOMIZABLE PROTOCOLS

Preloaded Protocols  
Gain access to the protocols the experts in Oral Immunotherapy use in their practices. With the OITConnect™ built-in protocols and workflows, you have the tools for the most efficient documentation for OIT treatment. The system is easy to use and modify for Oral Food Challenges.

OITCONNECT19:18:41Medical Clinic Inc.DANDAMUDI SESHAGIRI, Detroit (EST)

My ActivitiesPatient ChartMessagesReceptionScheduleSchedule ReportsFinalize NotesPractice Setup

Select PatientSearch Options

DashboardPractice Setup

GeneralBillingScheduleEMRPrivilegesRegister IP AddressOIT

SetupOIT

FoodsFood TypesSolutionsProtocolsReactionsMedicationsUnitsConcentrationsManufacturers

Search:Status: ActiveNew

NameActive

Almond (1)

AlmondYes

Cashew (1)

CashewYes

Cow's Milk (1)

Cow's Milk Horizon OrganicYes

Eggs (1)

Liquid Egg WhiteYes

Hazelnut (2)

Food ChallengeYes

HazelnutYes

Peanut (1)

PeanutYes

Sesame Seeds (1)

SesameYes

Sesame Tahini (1)

Sesame TahiniYes

Soy (1)

SoyYes

Walnut (1)

WalnutYes

Wheat (1)

WheatYes

Update Protocol

Name: Almond

Food: Almond

Active: ☒

Protocol Details - Almond

Seq#	Dosing	Location	Solutions / Food Types	Dose	Unit	Protein (mg)	Frequency	Duration(In Day)	Add	Delete
1	DAY ONE	In Office	Almond Solution A 1 mcg/ml	2.00	ml		every 20 minutes		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
2	DAY ONE	In Office	Almond Solution A 1 mcg/ml	4.00	ml		every 20 minutes		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
3	DAY ONE	In Office	Almond Solution B 10 mcg/ml	1.00	ml		every 20 minutes		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
4	DAY ONE	In Office	Almond Solution B 10 mcg/ml	2.00	ml		every 20 minutes		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
5	DAY ONE	In Office	Almond Solution B 10 mcg/ml	4.00	ml		every 20 minutes		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
6	DAY ONE	In Office	Almond Solution C 100 mcg/ml	1.00	ml		every 20 minutes		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
7	DAY ONE	In Office	Almond Solution C 100 mcg/ml	2.00	ml		every 20 minutes		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
8	DAY ONE	In Office	Almond Solution C 100 mcg/ml	4.00	ml		every 20 minutes		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
9	HOME DOSE	At Home	Almond Solution C 100 mcg/ml	4.00	ml		once a day	6	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
10	UPDOSE	In Office	Almond Milk	0.25	ml		once a day	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
11	HOME DOSE	At Home	Almond Milk	0.25	ml		once a day	6	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
12	UPDOSE	In Office	Almond Milk	0.50	ml		once a day	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
13	HOME DOSE	In Office	Almond Milk	0.50	ml		once a day	6	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
14	UPDOSE	In Office	Almond Milk	1.00	ml		once a day	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
15	HOME DOSE	At Home	Almond Milk	1.00	ml		once a day	6	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
16	UPDOSE	In Office	Almond Milk	1.50	ml		once a day	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
17	HOME DOSE	At Home	Almond Milk	1.50	ml		once a day	6	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
18	UPDOSE	In Office	Almond Milk	2.00	ml		once a day	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
19	HOME DOSE	At Home	Almond Milk	2.00	ml		once a day	6	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
20	UPDOSE	In Office	Almond Milk	3.00	ml		once a day	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
21	HOME DOSE	At Home	Almond Milk	3.00	ml		once a day	6	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
22	UPDOSE	In Office	Almond Milk	4.00	ml		once a day	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
23	HOME DOSE	At Home	Almond Milk	4.00	ml		once a day	6	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
24	UPDOSE	In Office	Almond Milk	6.00	ml		once a day	1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
25	HOME DOSE	At Home	Almond Milk	6.00	ml		once a day	6	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>



# UNITS and various concentrations – mcg, mg, grams, ml, oz, capsules

OITCONNECT Medical Clinic Inc., Detroit (EST)

Select Patient  Search Options

Dashboard Practice Setup X

General Billing Schedule EMR Privileges Register IP Address OIT

Setup OIT X

Foods Food Types Solutions Protocols Reactions Medications Units Concentrations Manufacturers

Search:  Status: Active  New

**Unit**

Unit	Active
<b>Dose (14)</b>	
5	Yes
mg	Yes
ml	Yes
tsp	Yes
tbsp	Yes
number	Yes
lb	Yes
microgram	Yes
grams	Yes
Slice	Yes
1	Yes
2	Yes
Nuts	Yes
Nut	Yes
<b>Concentrations (6)</b>	
mg/ml	Yes
mcg/ml	Yes
Stock	Yes
mcg	Yes
mg	Yes
grams/ml	Yes

OITCONNECT 19:21:44  
Medical Clinic Inc.  
DANDAMUDI SESHAGIRI, Detroit (EST)

Select Patient  Search Options

Dashboard Practice Setup X

General Billing Schedule EMR Privileges Register IP Address OIT

Setup OIT X

Foods Food Types Solutions Protocols Reactions Medications Units Concentrations Manufacturers

Search:  Status: Active  New

Name	Unit	Active
250	mcg/ml	Yes
25	mg/ml	Yes
2.5	mg/ml	Yes
0.25	mg/ml	Yes
0.025	mg/ml	Yes
0.0025	mg/ml	Yes
3	mg/ml	Yes
6	mcg/ml	Yes
12	mcg/ml	Yes
0.1	mg/ml	Yes
1	mg/ml	Yes
100	mcg/ml	Yes
Whole Milk	Stock	Yes
10	mg/ml	Yes
whole Chickpea	Stock	Yes
7	mg/ml	Yes
700	mcg/ml	Yes
70	mcg/ml	Yes
7	mcg/ml	Yes
200	Stock	Yes
5	mg/ml	Yes
Baked Milk Muffin	Stock	Yes
50	mg/ml	Yes
500	mcg	Yes
50	mcg	Yes
5	mcg	Yes
1.5	mg	Yes
250	mg/ml	Yes
25	mcg/ml	Yes
2.5	mcg/ml	Yes
0.25	mcg/ml	Yes
1.04	mg/ml	Yes
10	mcg/ml	Yes
1	mcg/ml	Yes

**Update Concentration**

Name: 250 mcg/ml

Active: ☒



# MANUFACTURERS – Customize for your own Inventory Management

Add in your own custom list of Manufacturers for better inventory control.

**OITCONNECT** 19:25:12  
Medical Clinic Inc.  
DANDAMUDI SESHAGIRI, Detroit (EST)

**My Activities**  
Patient Chart Messages Reception Schedule Schedule Reports Finalize Notes Practice Setup

Select Patient Search Options

**Dashboard Practice Setup**

General Billing Schedule EMR Privileges Register IP Address OIT

**Setup OIT**

Foods Food Types Solutions Protocols Reactions Medications Units Concentrations **Manufacturers**

Search: New Status: Active

Name	Address	State	City	Zip Code	Active
America's Peanut Farmers	8359 Office Park Drive	MI	Grand Blanc	48439	Yes
Byrd Mill	Testing Purpose	NY	Shirley	11967	Yes
Horizon Organic Vit D Whole Milk	From Grocery Store				Yes
Egglands' Best Liquid Egg Whites	Grocery Store				Yes
All Whites 100% Liquid Egg Whites	Grocery Store				Yes
Bush's Best Garbanzo Beans (Chic...	Grocery Store				Yes
Bob's Redmill Garbanzo Bean flour	Grocery store stone ground...				Yes
Bob's Red Mill Organic Coconut Fl...	Bobsredmill.com/organic-c...				Yes
Silk Organic Soy Milk (8gm protei...	Grocery Store				Yes
Sesame Seed (Unhulled)	Nuts.com				Yes
Sesame Tahini Spread (Max)	amazon.com				Yes

**Update Manufacturer**

Name: America's Peanut Farmers

Address: 8359 Office Park Drive

State: Michigan

City: Grand Blanc

ZIP Code: 48439

Active: ☒

# SCHEDULE



# OIT SCHEDULE

The schedule supports multiple locations, providers, and visit types. It is easy to schedule an appointment by double clicking on the appointment time.

**OITCONNECT** 20:05:20  
Medical Clinic Inc.  
DANDAMUDI SESHAGIRI, Detroit (EST)

**My Activities**  
Patient Chart Messages Reception Schedule Schedule Reports Finalize Notes Practice Setup

**OITConnect, Patient Te** Search Options

**Dashboard** **Schedule X**

View Office Visit Appointments View Calendar Multi Provider Day View View Facility Visit Appointments

**View Office Visit Appointments X**

Office Visit Appointments for - 4/12/2021 - Detroit - Giri Dandamudi, M.D. - All Visit Types (Please check the label)

Reminders Status Eligibility Patient Quick Search: Select Patient

Start Time	Visit Type	Duration (Min)	Auth- orization	Eligibility	Patien Portal	Added To No Show	Patient Name	Notes
8:30 AM	OIC							
8:45 AM	ODO							
9:00 AM	ODO							
9:15 AM	ODO	15	N		N		Mora, Michelle	
9:30 AM	UDO							
9:45 AM	UDO							
10:00 AM	UDO							
10:15 AM	ODO							
10:30 AM	ODO							
10:45 AM	UDO	15	N		N		New York, NJ	
11:00 AM	ORD							
11:15 AM	UDO							
11:30 AM	OGD							
1:00 PM	OIC	15	N		Y		Manor, Kris	
1:15 PM	UDO							
1:30 PM	UDO							
1:45 PM	UDO							
2:00 PM	UDO							
2:15 PM	UDO							
2:30 PM	UDO							
2:45 PM	MDO							

**New Appointment**

Existing Patient: Select Patient

Select Location: Detroit ☐ Wait List

Provider/User: Giri Dandamudi, M.D.

Referral Provider:

When: 4/12/2021 9:45 AM to 10:00 AM

Appointment Type: UPDOSE

Visit Reason:

Insurance Notes:

Special Instructions:

Ledger Comments:

Auth Info:

**Patient Details**

Account #:

Pt Balance:

Ins Balance:

Last Visit:

Insurance:

Next Appt:


**Earliest Appointments**

[New Patient Schedule...](#)

# PRE-LOADED and CUSTOMIZABLE SCHEDULE

As a general guide, we have preloaded Visit Types for your OIT schedule:

Initial Consult  
Day One  
Uprise  
Graduation Day  
Observation Dose  
Maintenance Dose  
Home Dose



Visit Type	Short Name	Instructions	Active	Default	Facility Visit	Tele Consult	Injection Visit
INITIAL CONSULT	OIC	Pre-vist template and check mark.	Yes	No	No	No	No
DAY ONE	ODO	Wait time 4-6 hours.	Yes	No	No	No	No
UPDOSE	UDO		Yes	No	No	No	No
GRADUATION DAY	OGD		Yes	No	No	No	No
OBSERVATION DOSE	OBD		Yes	No	No	No	No
MAINTENANCE DOSE	MDO		Yes	No	No	No	No
HOME DOSE	HDO		Yes	No	No	No	No



# SCHEDULE – PROVIDER SCHEDULE SETUP

Create custom OIT Schedules for your practice providers and locations. Here you can **Add Provider Schedule** and assign Visit Types to each of your providers preferred time slots.

**OITCONNECT** 19:31:54  
Medical Clinic Inc.  
DANDAMUDI SESHAGIRI, Detroit (EST)

My Activities  
Patient Chart Messages Reception Schedule Schedule Reports Finalize Notes Practice Setup

Select Patient Search Options Park Tabs Logout

Dashboard Practice Setup X

General Billing Schedule EMR Privileges Register IP Address OIT

Setup Schedule X

Visit Type Provider Schedule Setup Block Switch

Search

Search Criteria  
From Date: 06/15/2021  
To Date: 06/15/2021

Filter Criteria  
Select: ☒ Provider ☐ User  
Provider: Select Provider  
Location: Select Location  
Status: Active

Search Reset

Add Provider Schedule

Location	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday	Sunday
Schedule: 01/01/2021 - 12/31/2021 Provider: Clinician Lead Name:							
Telemedicine	08:00 AM - 08:45 PM	08:00 AM - 08:45 PM	08:00 AM - 08:45 PM	08:00 AM - 08:45 PM	08:00 AM - 08:45 PM	08:00 AM - 08:45 PM	08:00 AM - 08:45 PM
Edit Provider Schedule							
Schedule: 11/27/2020 - 11/27/2021 Provider: ATUL SHAH Name: New schedule Shirley							
Shirley	08:30 AM - 03:00 PM	08:30 AM - 03:00 PM	08:30 AM - 03:00 PM	08:30 AM - 03:00 PM	08:30 AM - 03:00 PM		
Edit Provider Schedule							
Schedule: 11/28/2020 - 11/28/2021 Provider: Douglas Jones Name: OIT schedule							
Layton	08:30 AM - 03:00 PM	08:30 AM - 03:00 PM	08:30 AM - 03:00 PM	08:30 AM - 03:00 PM	08:30 AM - 03:00 PM	08:30 AM - 03:00 PM	
Edit Provider Schedule							
Schedule: 01/04/2021 - 01/04/2022 Provider: Giri Dandamudi, M.D Name: DR.D							
Detroit	08:30 AM - 03:30 PM	08:30 AM - 03:30 PM	08:30 AM - 03:30 PM	08:30 AM - 03:30 PM	08:30 AM - 03:30 PM	08:30 AM - 03:30 PM	08:30 AM - 03:30 PM
Edit Provider Schedule							



# SCHEDULE – PROVIDER SCHEDULE SETUP

Select a **Provider**

**Name** the Schedule

Select the appointment **Durations**

Select a date range, **From Date & To Date**

Next, select the which days for each **Location**  
for which the schedule should be created.

Select **Next**.

Location	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday	Sunday
Detroit	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Telemedicine	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Dubai	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Select **Time** range for each day and location  
and click the **Green Arrow**.

Select the **Visit Type** for each time by clicking  
on the empty Visit Type field. To add a visit  
type to multiple times, click the checkbox for  
each time before selecting the Visit Type.

To Copy a schedule across days and locations,  
select the copy checkbox for the schedule you  
would like to copy.

When finished, select **Save**.

Time	Visit Type
6:00 AM	GRADUATION DAY
6:15 AM	GRADUATION DAY
6:30 AM	DAY ONE
6:45 AM	UPDOSE
7:00 AM	DAY ONE
7:15 AM	UPDOSE
7:30 AM	INITIAL CONSULT
7:45 AM	INITIAL CONSULT
8:00 AM	OBSERVATION DOSE
8:15 AM	INITIAL CONSULT
8:30 AM	MAINTENANCE DOSE
8:45 AM	GRADUATION DAY
9:00 AM	GRADUATION DAY
9:15 AM	GRADUATION DAY





# SCHEDULE – BLOCK AND SWITCH

Need to block a provider schedule for a vacation? Or maybe you need to switch appointments from one provider's schedule to another. With the Block and Switch setup, you can quickly do so.

The screenshot shows the OIT Connect Practice Setup interface for the 'Schedule' tab. The 'Block Setup' section is active, displaying a form for creating a schedule block. The form includes fields for 'Provider' (selected as 'Main Provider'), 'Location' (selected as 'Detroit'), 'From Date' (06/15/2021), 'To Date' (06/15/2021), and 'Block Name'. Below these fields is a table for 'Block Setup' with columns for 'Day', 'Start Time', and 'End Time'. The 'Day' column lists days from Monday to Sunday, each with a checkbox. The 'Start Time' and 'End Time' columns are empty. To the right of the 'Block Setup' section is the 'Appointment Setup Info' section, which shows 'Day Info' and a message 'No records found.'.

Day	Start Time	End Time
<input type="checkbox"/> Monday		
<input type="checkbox"/> Tuesday		
<input type="checkbox"/> Wednesday		
<input type="checkbox"/> Thursday		
<input type="checkbox"/> Friday		
<input type="checkbox"/> Saturday		
<input type="checkbox"/> Sunday		

Appointment Setup Info

Day Info

No records found.


The screenshot shows the OIT Connect Practice Setup interface for the 'Schedule' tab. The 'Switch' section is active, displaying a table for schedule switches. The table has columns for 'Schedule Date', 'From Provider', 'To Provider', 'Location', 'Switched By User', 'Switched Date', 'Undo By User', and 'Undo Date'. The table contains one record for the date 06/17/2021, where the schedule was switched from 'Main Provider' to 'Clinician Lead' at the 'Detroit' location, by 'Seshagiri Dandamudi'.

Schedule Date	From Provider	To Provider	Location	Switched By User	Switched Date	Undo By User	Undo Date
06/17/2021	Main Provider	Clinician Lead	Detroit	Seshagiri Dandamudi	06/15/2021		










# EMR – CARE MANAGEMENT



Create different categories and templates for Care Manager Notes so your practice can ensure documentation of care and consultations at every level.

 19:37:23

Medical Clinic Inc.  
DANDAMUDI SESHAGIRI, Detroit (EST)

My Activities


 Patient Chart  Messages  Reception  Schedule  Schedule Reports  Finalize Notes  Practice Setup

 Select Patient  Search Options

Dashboard Practice Setup X

General Billing Schedule EMR Privileges Register IP Address OIT

Care Management X

 New | Status: Active

Name	Characteristics	Discharge Date	Contact Date	Date Closed
General	Single Note			
Longitudinal care plan				
Episodic care plan				
Hospital Discharge Follow up	Series Note	Yes	Yes	Yes
ED Visit Followup	Series Note	Yes	Yes	Yes
Medication Understood	Single Note			
CPO	Series Note	Yes	Yes	



[www.oitconnect.com](http://www.oitconnect.com)



[www.oitconnect.com](http://www.oitconnect.com)



[www.oitconnect.com](http://www.oitconnect.com)

# OIT DASHBOARD



# DAILY PATIENT TRIAGE

The Daily Patient Triage tab is your main hub. View all patients scheduled, checked-in, and checked-out for the day.

You can check patient's in and access various areas of the patient's cart from here, including OIT Dose Administration.

19:44:17		<b>My Activities</b>													
Medical Clinic Inc. DANDAMUDI SESHAGIRI, Detroit (EST)		Patient Chart  Messages  Reception  Schedule  Schedule Reports  Finalize Notes  Practice Setup													
Select Patient  Search Options															
<b>Dashboard</b>															
<b>Daily Patient Triage</b>		<b>Office Communication</b> <b>Messages</b> <b>Patient Portal</b> <b>OIT Orders</b>													
Search															
<input type="checkbox"/>	Time	Check In Time	Visit Type	Room	Status	Provider	Order Id	Order Date	Order Status	Food	Image	Last Dose	Last Dose Time	Next Dose	Timer
<b>Patient: More, Michelle Account: 97775 DOB: 01/03/1984 (3)</b>															
<input type="checkbox"/>	09:15 AM	09:04 PM	ODO		Checked In	Giri Dandamudi, M.D	107	02/21/2021	Dose Admi...	Rye		Rye Flour Sol A 1 mg/m...	03/26/2021 10:54 AM	Rye Flour Sol A 1 mg/m...	
<input type="checkbox"/>	09:15 AM	09:04 PM	ODO		Checked In	Giri Dandamudi, M.D	79	01/22/2021	Dose Admi...	Peanut		PN Solution C 250 mcg/...	03/23/2021 1:14 PM	PN Solution C 250 mcg/...	
<input type="checkbox"/>	09:15 AM	09:04 PM	ODO		Checked In	Giri Dandamudi, M.D	73	01/04/2021	Dose Admi...	Cow's Milk		Cow's Milk Soln B 1 mg/...	03/23/2021 12:43 PM	Cow's Milk Soln B 1 mg/...	
<b>Patient: Manor, Kris Account: 97787 DOB: 03/03/1978 (1)</b>															
<input type="checkbox"/>	01:00 PM		OIC		Scheduled	Giri Dandamudi, M.D	113	03/03/2021	Dose Admi...	Cow's Milk		Cow's Milk Soln B 1 mg/...	03/08/2021 12:16 PM	Cow's Milk Soln B 1 mg/...	
<b>Patient: New york, NJ Account: 97799 DOB: 04/04/2016 (2)</b>															
<input type="checkbox"/>	10:45 AM		UDO		Scheduled	Giri Dandamudi, M.D	137	04/04/2021	Dose Admi...	Cow's Milk		Cow's Milk Soln B 1 mg/...	05/04/2021 1:37 PM	Cow's Milk Soln B 1 mg/...	
<input type="checkbox"/>	10:45 AM		UDO		Scheduled	Giri Dandamudi, M.D	138	04/04/2021	Dose Admi...	Peanut		PN Solution A 2.5 mcg/...	04/04/2021 9:14 AM	PN Solution B 25 mcg/...	



# OFFICE COMMUNICATION

Sometimes you need to communicate with other OITConnect users about non-patient related things. This Office Communication tool allows you to do just that.

Compose messages, attach files, and star important messages for later review.

The screenshot displays the OITConnect Office Communication interface. At the top, the header includes the OITConnect logo, the user's name (Medical Clinic Inc. DANDAMUDI SESHAGIRI, Detroit (EST)), the time (19:46:00), and a 'My Activities' section with links to Patient Chart, Messages, Reception, Schedule, Schedule Reports, Finalize Notes, and Practice Setup. Below the header is a 'Select Patient' dropdown and a 'Search Options' button. The main content area is titled 'Dashboard' and features tabs for 'Daily Patient Triage', 'Office Communication' (selected), 'Messages', 'Patient Portal', and 'OIT Orders'. Under the 'Office Communication' tab, there are sub-tabs for 'Inbox', 'Compose', 'Drafts', 'Sent', and 'Closed'. The 'Inbox' sub-tab is active, showing a list of messages. The list has columns for 'Select', 'Status', 'Priority', 'Date', 'From', 'To', and 'Subject'. A single message is listed with a status of 'Unread' (star icon), a priority of 'High' (exclamation mark icon), a date of '01/22/2021 04:11 PM', from 'DR.D Giri', to 'Admin Admin, Nurse Nurse...', and subject 'First Demo'. The interface also includes a 'Search Filter' on the left, a 'View Message' button on the right, and a footer showing 'Page 1 of 1' and 'Displaying 1 - 1 of 1'.


Select	Status	Priority	Date	From	To	Subject
<input type="checkbox"/>	☆	!	01/22/2021 04:11 PM	DR.D Giri	Admin Admin, Nurse Nurse...	First Demo



# MESSAGES

HIPAA Secure Messaging allows you to communicate with your patients through the OITConnect Patient Portal App.

Compose messages, read messages from patients, and forward messages to others OITConnect users. All Messages are attached to the patient's chart.



The screenshot displays the OITConnect Messages interface. At the top, there's a header with the OITConnect logo, user information (Medical Clinic Inc., DANDAMUDI SESHAGIRI, Detroit (EST)), and navigation tabs (Patient Chart, Messages, Reception, Schedule, Schedule Reports, Finalize Notes, Practice Setup). Below this is a search bar and a 'Select Patient' dropdown. The main content area is titled 'Inbox Messages' and contains a table of messages. The table has columns for Select, Status, Priority, Date, Account #, Patient Name, Patient Site, From, To, and Subject. The messages are listed in descending order of date. The bottom of the interface shows a pagination bar indicating 'Page 1 of 7' and a status bar at the very bottom showing 'Patient Portal' and 'Displaying 1 - 20 of 134'.

Select	Status	Priority	Date	Account #	Patient Name	Patient Site	From	To	Subject
<input type="checkbox"/>	☆	✉	3/21/2021 3:38 PM	97757	OITConnect, Pati...	Detroit	Patient	Doctor3 OITConne...	Requesting an Appointment
<input type="checkbox"/>	☆	✉	3/15/2021 2:48 AM	97757	OITConnect, Pati...	Detroit	Patient	Doctor3 OITConne...	Requesting an Appointment
<input type="checkbox"/>	☆	✉	3/6/2021 3:09 PM	97787	Manor, Kris	Detroit	Patient	Doctor3 OITConne...	Send vital information to practice
<input type="checkbox"/>	☆	✉	3/6/2021 9:07 AM	97787	Manor, Kris	Detroit	Patient	Doctor3 OITConne...	Send vital information to practice
<input type="checkbox"/>	☆	✉	3/5/2021 6:28 PM	97787	Manor, Kris	Detroit	Patient	Doctor3 OITConne...	Requesting an Appointment
<input type="checkbox"/>	☆	✉	3/5/2021 6:28 PM	97787	Manor, Kris	Detroit	Patient	Doctor3 OITConne...	Send vital information to practice
<input type="checkbox"/>	☆	✉	3/5/2021 4:45 PM	97787	Manor, Kris	Detroit	Patient	Doctor3 OITConne...	Send vital information to practice
<input type="checkbox"/>	☆	✉	3/5/2021 10:34 AM	97757	OITConnect, Pati...	Detroit	Patient	Doctor3 OITConne...	Requesting an Appointment
<input type="checkbox"/>	☆	✉	3/5/2021 8:08 AM	97788	Gibbs, Jack	Detroit	Patient	Doctor3 OITConne...	Successful access to patient portal
<input type="checkbox"/>	☆	✉	3/4/2021 6:48 PM	97787	Manor, Kris	Detroit	Patient	Doctor3 OITConne...	
<input type="checkbox"/>	☆	✉	3/3/2021 8:39 PM	97787	Manor, Kris	Detroit	Patient	Doctor3 OITConne...	Need appointment
<input type="checkbox"/>	☆	✉	3/3/2021 8:25 PM	97787	Manor, Kris	Detroit	Patient	Doctor3 OITConne...	Successful access to patient portal
<input type="checkbox"/>	☆	✉	3/3/2021 6:05 AM	97778	Force, Andrea	Detroit	Patient	Doctor3 OITConne...	Successful access to patient portal
<input type="checkbox"/>	☆	✉	2/22/2021 11:29 AM	97757	OITConnect, Pati...	Detroit	Patient	Doctor3 OITConne...	Doctor
<input type="checkbox"/>	☆	✉	2/19/2021 2:32 PM	97757	OITConnect, Pati...	Detroit	Patient	Doctor3 OITConne...	next appointment
<input type="checkbox"/>	☆	✉	2/19/2021 3:12 AM	97757	OITConnect, Pati...	Detroit	Patient	Doctor3 OITConne...	demo
<input type="checkbox"/>	☆	✉	2/19/2021 3:11 AM	97757	OITConnect, Pati...	Detroit	Patient	Doctor3 OITConne...	demo
<input type="checkbox"/>	☆	✉	2/19/2021 3:11 AM	97757	OITConnect, Pati...	Detroit	Patient	Doctor3 OITConne...	sample
<input type="checkbox"/>	☆	✉	2/19/2021 2:35 AM	97757	OITConnect, Pati...	Detroit	Patient	Doctor3 OITConne...	others
<input type="checkbox"/>	☆	✉	2/19/2021 2:34 AM	97757	OITConnect, Pati...	Detroit	Patient	Doctor3 OITConne...	for patient



# PATIENT PORTAL

Sign-up patients for the OITConnect Patient Portal so they can check-in for Office Visits, document their at-home doses, and communicate directly with your practice. It is simple to register and manage patients accounts from a centralized location.

**OITCONNECT** 19:47:29  
Medical Clinic Inc.  
DANDAMUDI SESHAGIRI, Detroit (EST)

My Activities  
Patient Chart Messages Reception Schedule Schedule Reports Finalize Notes Practice Setup

Select Patient Search Options

Dashboard  
Daily Patient Triage Office Communication Messages Patient Portal OIT Orders

Register Web User

Account #	Name	Date of Birth	Cell Phone	Registered Date	User Type	Status	Login	Delete
97786	Jimmy Carter	03/03/1980		03/03/2021	Authorized User	✓	🔒	🗑️
User: Kris Manor Email: giridandam@gmail.com								
97787	Kris Manor	03/03/1978	2699683030	03/03/2021	Patient	✓	🔒	🗑️
User: DemoOffice Manager Email: mgr1.oitconnect@gmail.com								
1463	Patient Ann 4	06/12/1969		09/11/2020	Patient	✓	🔒	🗑️
User: Patient 106 Email: nouser1@outlook.com								
73845	Patient D 106	02/03/2007	111-212-1212	11/13/2020	Patient	✓	🔒	🗑️
75810	Patient H 114	05/27/2009		11/13/2020	Patient	✓	🔒	🗑️
User: Patient 102 Email: nouser@outlook.com								
71422	Patient George 102	02/18/1977	810-650-1401	11/13/2020	Patient	✓	🔒	🗑️
User: Twenty March Email: olivehealthglobal.com								
97793	Twenty March	03/13/2014		03/13/2021	Patient	✓	🔒	🗑️
User: Jack Gibbs Email: olivehealthglobal@gmail.com								
97788	Jack Gibbs	03/05/1999	2699683030	03/03/2021	Patient	✓	🔒	🗑️
User: Patient OITConnect Email: patient1.oitconnect@gmail.com								
97757	Patient Test OITConnect	07/13/1972	088-424-5239	09/11/2020	Patient	✓	🔒	🗑️
User: Jane Doe Email: patient1.sipproject@gmail.com								
71404	Jane Doe	06/01/1999		11/13/2020	Patient	✓	🔒	🗑️
User: Patient Health Email: patientrevival2@gmail.com								
97789	Patient Patient Health			03/03/2021	Patient	✓	🔒	🗑️

Page 1 of 1

**Register Web User**

Save Print Consent Form

Patient:\* Select patient Relation:\* Self

Name:\* First Name Last Name Suffix

Address: Address 1 Address 2

City, State, ZIP: City Select State Zip

Country: USA

Phone: Home Phone Cell Phone

Email:\* Email

# OIT WORKFLOW AND HOME DOSES

- **90 % of OIT is done at home**
- **We have that covered!**


# DEMOGRAPHICS

Capture your patient demographic information. With an EMR integration, this information will pre-populate, so there will be no need for duplicate data entry.

Check out our integrated EMR partner

**ModuleMD**<sup>®</sup>  
Providing solutions...not just software

or ask about an integration with your existing system.



The screenshot displays the OITCONNECT Patient Demographics form. At the top, the header includes the OITCONNECT logo, the time 19:59:02, and the user's name Medical Clinic Inc. DANDAMUDI SESHAGIRI, Detroit (EST). Below the header, there are tabs for Patient Chart, Messages, Reception, Schedule, Schedule Reports, Finalize Notes, and Practice Setup. The main form area is titled "Patient Demographics" and contains several sections: "Patient Information" with fields for Last Name, First Name, Middle Name, Suffix, Nick Name, Birth Sex, Date of Birth, Race, Ethnicity, Preferred Language, Primary Language, Occupation, Marital Status, Account Status, Updated On, Statement Preference, SSN, Employer, and Preliminary Cause of Death; "Patient Communication Information" with fields for Suite/Apt, City, State, ZIP Code, Email, Cell Phone, Home Phone, Work Phone, Ext, Communication Preference, Number Of Reminders, Emergency Contact, and Relation to Patient; and "Patient Provider Information and Pharmacy Information" with fields for Location, Facility, Primary Care Provider, Referral Provider 1, Referral Provider 2, Referral Source, Pharmacy Retail, Pharmacy Mail Order, Address, Address2, City, State, and Zip.



# MEDICATIONS

Medications such as EMT/ER Intervention, Inhalers, and Oral Medications can be added to the system repository.

**OITCONNECT** 19:20:04  
Medical Clinic Inc.  
DANDAMUDI SESHAGIRI, Detroit (EST)

My Activities  
Patient Chart Messages Reception Schedule Schedule Reports Finalize Notes Practice Setup

Select Patient Search Options

Dashboard Practice Setup X

General Billing Schedule EMR Privileges Register IP Address OIT

Setup OIT X

Foods Food Types Solutions Protocols Reactions Medications Units Concentrations Manufacturers

Search: Status: Select All New

Medication	Strength	Active
Additional EMT/ER Intervention (1)		
Epinephrine Brand	43.00 mg	Yes
Inhaler or Nebulizer Medication (1)		
Inhaler-bronchodilator	10.00 5	Yes
Oral Medications (1)		
Antihistamine Brand	40.00 mg	Yes

Update Medication

Medication: Epinephrine Brand

Strength: 43.00 mg

Group: Additional EMT/ER Inte

Active: ☒

Update Clear Close

Dashboard Practice Setup X

General Billing Schedule EMR Privileges Register IP Address OIT

Setup OIT X

Foods Food Types Solutions Protocols Reactions Medications Units Concentrations Manufacturers

Search: Status: Select All New

Medication	Strength	Active
Additional EMT/ER Intervention (1)		
Epinephrine Brand	43.00 mg	Yes
Inhaler or Nebulizer Medication (1)		
Inhaler-bronchodilator	10.00 5	Yes
Oral Medications (1)		
Antihistamine Brand	40.00 mg	Yes

Update Medication

Medication: Epinephrine Brand

Strength: 43.00 mg

Group: Additional EMT/ER Inte

Active: ☒


# DOSE ADMINISTRATION: DAY 1 PATIENT TRIAGE

To administer a patient dose once a patient has been checked in for their appointment you must first confirm the order:

## Dashboard > Daily Patient Triage

Select the **checkbox** next to the order.

Select **Confirm Order**.




Time	Check In Time	Visit Type	Room	Status	Provider	Order Id	Order Date	Order Status	Food	Image	Last Dose	Last Dose Time	Next Dose	Timer
10:30 AM	09:09 AM	ODO		Checked Out	Giri Dandamudi, M.D	104	02/20/2021	Dose Admin...	Chickpea		Chickpea Sol A 10 mcg/...	03/27/2021 10:50 AM	Chickpea Sol A 10 mcg/...	
10:30 AM	09:09 AM	ODO		Checked Out	Giri Dandamudi, M.D	82	01/22/2021	Dose Admin...	Chickpea		Chickpea Fragment who...	6:14 AM	Chickpea whole (1.00)	
10:30 AM	09:09 AM	ODO		Checked Out	Giri Dandamudi, M.D	102	02/20/2021	Ordered	Cow's Milk				Cow's Milk Soln A 100 ...	
10:30 AM	09:09 AM	ODO		Checked Out	Giri Dandamudi, M.D	103	02/20/2021	Ordered	Eggs				Liq Egg White Solution ...	
10:30 AM	09:09 AM	ODO		Checked Out	Giri Dandamudi, M.D	146	06/15/2021	Ordered	Peanut				PN Solution A 2.5 mcg/...	

You will be asked to confirm the order. At this time, you can:

- Enter **Height & Weight**
- Check if the patient has a **Signed Consent Form** and has completed **Pre-Challenge Preparation**
- Make changes to the schedule

Once you have completed confirming the order, select **Update**.



Seq#	Dosing	Location	Solutions / Food Types	Dose	Unit	Protein (mg)	Frequency	Duration(In Days)	Add	Delete
1	DAY ONE	In Office	Cow's Milk Soln A 100 mcg/ml	2.00	ml	0.20	every 20 minutes	1		
2	DAY ONE	In Office	Cow's Milk Soln A 100 mcg/ml	4.00	ml	0.40	every 20 minutes	1		
3	DAY ONE	In Office	Cow's Milk Soln A 100 mcg/ml	6.00	ml	0.60	every 20 minutes	1		
4	DAY ONE	In Office	Cow's Milk Soln B 1 mg/ml	0.60	ml	0.60	every 20 minutes	1		
5	DAY ONE	In Office	Cow's Milk Soln B 1 mg/ml	0.80	ml	0.80	every 20 minutes	1		
6	DAY ONE	In Office	Cow's Milk Soln B 1 mg/ml	1.00	ml	1.00	every 20 minutes	1		
7	DAY ONE	In Office	Cow's Milk Soln B 1 mg/ml	2.00	ml	2.00	every 20 minutes	1		
8	DAY ONE	In Office	Cow's Milk Soln B 1 mg/ml	4.00	ml	4.00	every 20 minutes	1		
9	DAY ONE	In Office	Cow's Milk Soln B 1 mg/ml	6.00	ml	6.00	every 20 minutes	1		
10	HOME DOSE	At Home	Cow's Milk Soln B 1 mg/ml	6.00	ml	6.00	once a day	6		



# INTEGRATED REACTIONS PAGE WITH PRACTICE AND PATIENT APP

**OIT CONNECT** 19:19:28  
Medical Clinic Inc.  
DANDAMUDI SESHAGIRI, Detroit (EST)

**My Activities**  
Patient Chart Messages Reception Schedule Schedule Reports Finalize Notes Practice Setup

Select Patient Search Options Park Tabs Logout

**Dashboard Practice Setup X**

General Billing Schedule EMR Privileges Register IP Address OIT

**Setup OIT X**

Foods Food Types Solutions Protocols Reactions Medications Units Concentrations Manufacturers

Search: Status: Active New

Symptoms	Severity	Active
<b>Gut (2)</b>		
Mild nausea or discomfort	Mild	Yes
Repetitive vomiting	Severe	Yes
<b>Heart (4)</b>		
Dizziness	Severe	Yes
Weak pulse	Severe	Yes
Faintness	Severe	Yes
Pale or bluish skin	Severe	Yes
<b>Lung (3)</b>		
Repetitive cough	Severe	Yes
Wheezing	Severe	Yes
Shortness of breath	Severe	Yes
<b>Mouth (1)</b>		
Significant swelling of the tongue or lips	Moderate	Yes
<b>Nose (1)</b>		
Itchy or runny nose	Mild	Yes
<b>Skin (2)</b>		
A few hives	Mild	Yes
Many hives over body	Severe	Yes
<b>Throat (2)</b>		
Trouble breathing or swallowing	Severe	Yes
Tight or hoarse throat	Severe	Yes

**Update Reaction Template** Update Clear Close

Symptoms: Mild nausea or discomfort  
Severity: Mild  
Group: Gut  
Active: ☒

Symptoms	Severity	Active
<b>Gut (2)</b>		
Mild nausea or discomfort	Mild	Yes
Repetitive vomiting	Severe	Yes
<b>Heart (4)</b>		
Dizziness	Severe	Yes
Weak pulse	Severe	Yes
Faintness	Severe	Yes
Pale or bluish skin	Severe	Yes
<b>Lung (3)</b>		
Repetitive cough	Severe	Yes
Wheezing	Severe	Yes
Shortness of breath	Severe	Yes
<b>Mouth (1)</b>		
Significant swelling of the tongue or lips	Moderate	Yes
<b>Nose (1)</b>		
Itchy or runny nose	Mild	Yes
<b>Skin (2)</b>		
A few hives	Mild	Yes
Many hives over body	Severe	Yes
<b>Throat (2)</b>		
Trouble breathing or swallowing	Severe	Yes
Tight or hoarse throat	Severe	Yes



# OIT SOAP NOTE

It is easy to document in the OIT SOAP Note. With pre-built fields and templates for OIT consultations you can quickly watch your note build.

**OITCONNECT** 19:54:24  
Medical Clinic Inc.  
DANDAMUDI SESHAGIRI, Detroit (EST)

My Activities  
Patient Chart Messages Reception Schedule Schedule Reports Finalize Notes Practice Setup

Search Options Patient: OITConnect, Patient Test Account #: 97757 DOB: 07/13/1972 Age: 48 yrs, 11 mo Sex: M Home

Dashboard Patient Chart X

Demographics View Notes Visit Notes Messages OIT

Patient Dashboard OIT SOAP Notes X Change

To Do S O A P MU All

All Categories Deselect All

- Chief Complaint
- Physical Examination
- Vital Signs
- Systemic Examination
- Assessment
- Plan
- Patient Instructions

Visit Slip: 961648 Visit Date: 06/07/2021 Provider: Giri Dandamudi, M.D.

Memo Encounter Form Copy History Print Fax Flag To Delete PIN Approve << Previous Next >> Review CDS Visit Summary

**Chief Complaint**  
Patient is in today for initial consult for peanut allergy. Had 3 systemic reactions. Had 2 ER visits. He wants to start OIT program. New template

**Physical Examination**  
Vital signs: Height: 5.00 feet 10 inches. Weight: 189.0 pounds. BMI: 27.1. BP: Sitting 115/78 mmHg. Respiratory rate: 15/min. PEFR: 95 L/min. Pulse O<sub>2</sub>: 99.0 %. ACT/ATAQ: 21. Inhaled Oxygen Concentration: 99.0 %. FeNO: 37 ppb. UCT: 5.  
Systemic Examination: EYES: Conjunctivae clear bilaterally, with no injection or edema. There is no discharge. Lids are normal with no swelling or redness.  
EARS: Canals are clear with no cerumen. Tympanic membranes are transparent with no fluid or erythema present.  
NOSE: Turbinates are pink with no swelling. There is no discharge present. There are no polyps, masses.  
PHARYNX: Pharynx is normal with no exudates or postnasal drip.  
LUNGS: No signs of respiratory distress, no retractions or accessory muscle use.  
LUNGS: Lungs are clear to auscultation bilaterally. There are no rhonchi, wheezing, or rales.  
HEART: Heart has regular rate and rhythm with no murmur.  
SKIN: No hives or eczema

**Assessment**

ICD Code	SNOMED	Code Description	Created Date	Comments
Z91.02		Food additives allergy status	6/15/2021	
Z91.018		Allergy to other foods	6/15/2021	

**Plan**  
1. Initial consult for peanut is done.  
2. All consent forms signed.  
3. Discussed benefits and risks with OIT peanut.  
4. Agree with financial billing and payments.  
5. Agreed to sign in for patient portal.  
6. Can start OIT anytime and I will follow up in 3 months.

**Patient Instructions**  
Keep contact number of provider all the time.  
Use patient portal and document all information. Look for notifications.  
Return visit as schedule and while in parking lot check in and wait.

**Chief Complaint**

Template Name: Search Kit

Initial Consult: IC

New temp: NT

Patient is in today for initial consult for peanut allergy. Had 3 systemic reactions. Had 2 ER visits. He wants to start OIT program. New template

Save Save As Template

**Vital Signs**

Height: 5 10 (Feet) (Inches) Weight: 189 (Pounds) (Ounces) Gms Kgs Lbs BMI: 27.1

Pulse rate: /min Rhythm: Select One Temperature: °C °F

BP sitting: Systolic 115 /mmHg Diastolic 78 /mmHg

BP standing: Systolic /mmHg Diastolic /mmHg

BP lying: Systolic /mmHg Diastolic /mmHg

Respiratory rate: 15 /min PEFR: 95 L/min Personal Best PEFR: L/min Pulse O<sub>2</sub>: 99.0 % ACT/ATAQ: 21 FeNO: 37 ppb UCT: 5

Head circumference: Cms Inhaled Oxygen Concentration: 99.0 %

General appearance:





# CARE MANAGER NOTES

You can set up templates for Care Manager Notes. Document everything from Contact Date, who Initiated the call, Duration, Facility, to Care Type.

**Patient Chart > Visit Notes > Care Manager Notes**

The screenshot displays the OITCONNECT Patient Chart interface for a patient named Demo, Mark. The interface includes a patient profile, a list of notes, and a 'New' button to add a note. The 'New' form includes fields for Date of Service, Note Type, Template, Notes, Discharge Date, Contact Date, Initiative, Duration, Facility, and Care Type.

**General (Single Note)**

Note Date	Location	Created By	Initiative By	Created Date	Updated Date	Discharge Date	Contact Date	Date Closed	Duration (Mins)	Facility	Care Type	Series Count
06/15/2021	Detroit	Seshagiri Dandamudi	Practice	06/15/2021					20			

**Add New**

Finalize Save As Draft Clear

Date of Service: 06/15/2021 Note Type: General (Single Note)

Template: Select

Notes: Called patient to consult

Discharge Date: Select Date... Contact Date: Select Date...


Initiative: Practice Duration: 20 Mins

Facility: Select Facility Care Type: Select Care Type

# DOSE ADMINISTRATION: DAILY PATIENT TRIAGE

To administer a dose, return to the Daily Patient Triage Dashboard. **Dashboard > Daily Patient Triage**


Select the **checkbox** next to the order.  
Select **Dose Administer**.




The screenshot shows the OITCONNECT interface for a patient named Demo, Mark. The 'Daily Patient Triage' tab is active. A table lists several orders, with the last one (Order ID 102) selected. The 'Dose Administer' button is highlighted with a yellow box.

Time	Check In Time	Visit Type	Room	Status	Provider	Order Id	Order Date	Order Status	Food	Image	Last Dose	Last Dose Time	Next Dose	Timer
<input type="checkbox"/>	10:30 AM	09:09 AM	ODO	Checked Out	Giri Dandamudi, M.D	104	02/20/2021	Dose Admi...	Chickpea		Chickpea Sol A 10 mcg/...	03/27/2021 10:50 AM	Chickpea Sol A 10 mcg/...	
<input type="checkbox"/>	10:30 AM	09:09 AM	ODO	Checked Out	Giri Dandamudi, M.D	82	01/22/2021	Dose Admi...	Chickpea		Chickpea Fragment who...	6:14 AM	Chickpea whole (1.00)	
<input type="checkbox"/>	10:30 AM	09:09 AM	ODO	Checked Out	Giri Dandamudi, M.D	103	02/20/2021	Ordered	Eggs				Liq Egg White Solution ...	
<input type="checkbox"/>	10:30 AM	09:09 AM	ODO	Checked Out	Giri Dandamudi, M.D	146	06/15/2021	Ordered	Peanut				PN Solution A 2.5 mcg/...	
<input checked="" type="checkbox"/>	10:30 AM	09:09 AM	ODO	Checked Out	Giri Dandamudi, M.D	102	02/20/2021	Confirmed	Cow's Milk				Cow's Milk Soln A 100 ...	

To administer a dose:

Select the **checkbox** under **Double Check**  
Select the **Green Arrow**  to Administer.



The screenshot shows the OITCONNECT 'Dose Administer' form for the selected order. The 'Double Check' checkbox and the 'Green Arrow' button are highlighted with yellow boxes.

**Account #:** 97757  
**Name:** OITConnect, Patient Test  
**DOB:** 07/13/1972  
**Sex:** Male

**Communication Preferences:** Email: patient1.oitconnect@gmail.com  
**Emergency Person:** last name, emergency contact  
**Relationship:**  
**Emergency Contact:**

**Provider:** Dandamudi, Giri  
**Cell Phone:**  
**Work Phone:**  
**Email:** patient1.oitconnect@gmail.com

**Administered By:** Dandamudi, Seshagiri  
**Administered Date:** 06/15/2021 8:25 PM  
**Consent Form Signed:** Yes  
**Pre Challenge Preparation Done:** Yes


**Check In Time:** 8:25 PM  
**Start Time:** 8:25 PM  
**Height:** 5 Feet 10 Inches  
**Weight:** 189 Lbs

Seq#	Dosing	Location	Solutions / Food Types	Dose	Unit	Protein (mg)	Dose Given	Date	Time	Double Check	Vitals	Reaction	Comments	
1	DAY ONE	In Office	Cow's Milk Soln A 100 mcg/ml	2.00	ml	0.20	2.00	06/15/2021	8:27 PM	<input checked="" type="checkbox"/>				

# INTEGRATED PATIENT APP



# PATIENT PORTAL - INTEGRATED



From the patient app, under **Oral Immunotherapy** in the menu, Patients can document at-home doses including **vitals**, **reactions**, and **comments**. Once they hit submit, the does administration will be documented and visible in the Provider's trending report.





# PATIENT PORTAL – PRACTICE PORTAL - INTEGRATED



Download on the  
**App Store**





**Google Play**



## Appointment Management

The OITConnect™ patient application tracks previous and **upcoming in-office dosing appointments**.

Patients can also schedule their next appointment with their provider through the application.







## Continued Communication

Patients can stay connected to the practice from home with a **two-way secure messaging system**. Engaged patients experience more effective OIT treatment. Providers and patients can monitor progress and ask questions through the OITConnect™ application.

# INTEGRATED PRACTICE AND PATIENT PORTAL



 **Dosing Reminders**  
Patients can stay on track with **at-home dosing through the OITConnect™ patient application.** Patients will be reminded of at-home maintenance dosing schedules.






# INTEGRATED PRACTICE AND PATIENT PORTAL



## Reaction Documentation

The OITConnect™ patient application makes it easier to **track patient progress during at-home maintenance dosing**. With standardized, easy to understand questionnaires, patients can report any minor reactions or let their providers know things are going well. This information is automatically sent to the provider, attached to the patient's chart, and visible in the patient OIT treatment trending.





# Trending Report

To view a trending report of doses administered, including at home doses documented in the patient app:

**Patient Chart > OIT > Orders**

Select the **checkbox** next to the order.  
Select **View**.

Order Id	Order Date	Appointment Date	Order Status	Account #	Patient Name	Date Of Birth	Food	Image	Provider	Location	Confirmed By	Confirmed Date	Finalized By	Finalized Date
146	06/15/2021	06/07/2021	Ordered	97757	OITConnect, Patient	07/13/1972	Peanut		Giri Dandamudi, M.D	Detroit				
<input checked="" type="checkbox"/> 145	06/03/2021	06/07/2021	Dose Administered	97757	OITConnect, Patient	07/13/1972	Peanut		Giri Dandamudi, M.D	Detroit	Dandamudi, Seshagiri	06/03/2021 11:40 AM		

**Order#:** 145 **Order Date:** 06/03/2021 **Location:** Detroit

**Patient**  
Account#: 97757  
Name: OITConnect, PatientTest  
DOB: 07/13/1972  
Sex: Male  
Order Confirmed  
Confirmed By: Dandamudi, Seshagiri  
Confirmed Date: 06/03/2021  
Consent Form Signed: No  
Pre Challenge Preparation Done: No

**Communication Preferences:** Email  
Email : patient1.oitconnect@gmail.com  
Emergency Person: last name, emergency contactmiddleham Relationship:  
Emergency Contact:

**Check In Time:** 11:40 AM  
**Start Time:** 11:40 AM  
**Height:**  
**Weight:**

**Provider**  
Provider: Dandamudi, GiriGiri  
Cell Phone:  
Work Phone:  
Email: giridandam@gmail.com  
**Order Finalized**  
Observation Time:  
End Time:  
Post Challenge Instructions: No  
Last Comments of the Day:

Seq	Dosing	Location	Solutions / Food Types	Dose	Protein (mg)	Frequency	Duration(In Days)	Dose Given	Given Time	Temperature	Pulse Rate	BP Systolic	BP Diastolic	R.R	Pulse Ox	FEV1
1	DAY ONE	In Office	PN Solution A 2.5 mcg/ml	2.00 ml		every 20 minutes	0	2.00	06/03/2021 11:40 AM							
2	DAY ONE	In Office	PN Solution A 2.5 mcg/ml	4.00 ml		every 20 minutes	0	4.00	06/03/2021 11:41 AM							
3	DAY ONE	In Office	PN Solution B 25 mcg/ml	1.00 ml		every 20 minutes	0	1.00	06/03/2021 11:41 AM							
4	DAY ONE	In Office	PN Solution B 25 mcg/ml	2.00 ml		every 20 minutes	0	2.00	06/03/2021 11:41 AM							
5	DAY ONE	In Office	PN Solution B 25 mcg/ml	4.00 ml		every 20 minutes	0	4.00	06/03/2021 11:41 AM							
6	DAY ONE	In Office	PN Solution C 250 mcg/ml	1.00 ml		every 20 minutes	0	1.00	06/03/2021 11:41 AM							
7	DAY ONE	In Office	PN Solution C 250 mcg/ml	2.00 ml		every 20 minutes	0	2.00	06/03/2021 11:41 AM							
8	DAY ONE	In Office	PN Solution C 250 mcg/ml	4.00 ml		every 20 minutes	0	4.00	06/03/2021 11:41 AM							
9	HOME DOSE	At Home	PN Solution C 250 mcg/ml	4.00 ml		once a day	6	4.00	06/03/2021 11:45 AM	100.00	98	100	80	20	98	400.00
10	HOME DOSE	At Home	PN Solution C 250 mcg/ml	4.00 ml		once a day	6	4.00	06/03/2021 11:50 AM	89.00	90	100	80			
11	HOME DOSE	At Home	PN Solution C 250 mcg/ml	4.00 ml		once a day	6	4.00	06/04/2021 11:29 AM							
12	HOME DOSE	At Home	PN Solution C 250 mcg/ml	4.00 ml		once a day	6	4.00	06/06/2021 07:46 PM	97.00	86	110	80	14	98	400.00
13	HOME DOSE	At Home	PN Solution C 250 mcg/ml	4.00 ml		once a day	6	4.00	06/15/2021 08:52 PM	96.00	80	100	80		95	470.00



Email:  
[Support@OITConnect.com](mailto:Support@OITConnect.com)

Sign Up:  
[www.OITConnect.com](http://www.OITConnect.com)



# OIT CONNECT

A G F T C O M P A N Y

Douglas H. Jones, MD, FACAAI, FAAAAI  
Atul N. Shah, MD, FACAAI, FAAAAI  
Giri Dandamudi, MD, FACAAI, FAAAAI  
Co-Founders

[www.oitconnect.com](http://www.oitconnect.com)